



Catalase Fluorometric Detection Kit

Product Information

Cat.No.

Kit-0180

Product Overview

A sensitive assay that utilizes a non-fluorescent detection reagent to measure H₂O₂ substrate as a result of the reaction of catalase. The production of hydrogen peroxide in eukaryotic cells is an end product of various oxidases and superoxide dismutase reactions. Accumulation of H₂O₂ can result in cellular damage through oxidation of proteins, DNA and lipids thus resulting in cell death and mutagenesis.

Applications

Fluorescence plate reader

Usage

1. For Research use only. Not for use in diagnostic procedures. 2. Practice safe laboratory procedures by wearing protective clothing and eyewear. 3. The fluorescent product of the detection reagent is not stable in the presence of thiols (DTT or 2-mercaptoethanol). Keep these reactants below 5mM.

Storage

1. The kit contains multiple storage temperature components. Please see labels of individual components for storage instructions. 2. Once a vial of the Detection reagent is opened, it should be used promptly and frozen since it is subject to oxidation by air.

Kit Components

Reagent-Storage Temperature 1. 5X Reaction Buffer pH 7.4: 25mL, 2-8°C; 2. Detection Reagent: 1 Vial, 2-8°C-Undiluted, Diluted: -20°C (Aliquot in Single Use Vials); 3. Horseradish Peroxidase: 1 Vial, 2-8°C-Undiluted, Diluted: -20°C; 4. Hydrogen Peroxide: 1 Vial, 2-8°C; 5. Catalase Enzyme: 1 Vial, 2-8°C (Crystalline Suspension)

Features & Benefits



CREATIVE **BIOMART**[®]
Assay Kit

Catalase Fluorometric Detection Kit

1. Sensitive Fluorescent Assay.2. Enzyme Positive Control included in kit.3. Can monitor multiple time points to follow kinetics.4. One-step, No Wash assay.5. Adaptable for High Throughput format.6. Monitors enzymatic activity.7. Applications-Fluorescent Plate Reader.

Tel: 1-631-559-9269 1-516-512-3133

Email:info@creative-biomart.com

Fax:1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA