

Fructose-6-Phosphate Assay Kit

Product Information

Cat.No.

Kit-0331

Product Overview

The PicoProbe F6P Assay Kit can detect F6P in the range of 0.01 to 0.5 nmoles with detection sensitivity ~1 μ M of F6P. (100 assays) Simple procedure; takes ~ 35-40 minutes to complete the assayFast and convenientSensitive, efficient and accurate assay

Size

100 assays

Applications

The PicoProbe F6P Assay Kit can detect F6P in the range of 0.01 to 0.5 nmoles with detection sensitivity ~1 μ M of F6P.

Target Species

Mammalian

Usage

This product is furnished for LABORATORY RESEARCH USE ONLY. Not for diagnostic or therapeutic use.

Kit Components

F6P Assay BufferF6P ProbeF6P Enzyme MixF6P ConverterF6P DeveloperF6PStandard (10 μ mol)

Detection method Fluorescence Plate ReaderFluorometer

Compatible Sample Types

Cell Culture SupernatantsFood ProductsGrowth MediaPlasmaSerumUrine