



Cell Invasion Assay Kit (Laminin), 8 ?m

Product Information

Cat

Kit-1001

Common Name

Cell

Cat.No.

Kit-1001

Description

Cell invasion is the ability of cells to migrate from one area to another through an extracellular matrix. Cell invasion is exhibited by both normal cells as well as cancerous cells in response to specific external signals, including chemical & mechanical stimuli. During invasion, extracellular matrix is enzymatically degraded by cellular proteases before cells migrate to the new location. Cell invasion is required for normal processes such as wound repair, vasculature formation and the inflammatory response as well as the abnormal invasion of tissues by tumor cells during metastasis. Cell Invasion Assay Kit utilizes a Boyden chamber coated with Laminin, where the cells invade the matrix and then migrate through a semipermeable membrane in the Boyden chamber in response to stimulants or inhibitory compounds. The percent cell invasion can be analyzed directly in a plate reader. Our assay is easy to use, sensitive and adaptable to high-throughput systems.

Applications

Measure cell invasion in response to stimuli.
Screen, study, or characterize compounds that influence cell invasion.

Storage

-20°C

Shipping

Gel Pack



Cell Invasion Assay Kit (Laminin), 8 ?m

Size

24 assays

Kit Components

Wash Buffer; Cell Dissociation Solution; Control Invasion Inducer; Cell Dye; Cell Invasion Chamber; Laminin

Detection method Fluorescence (Ex/Em = 530/590 nm)

Features & Benefits

Highly sensitive fluorometric method to measure cell invasion in response to a variety of biochemical stimuli;

Simple & High throughput-adaptable;

Reproducible, Quantitative tool for screening, studying, and characterizing compounds that affect cell invasion
