



# TNKS1 Histone Ribosylation Assay Kit (Biotin-labeled NAD<sup>+</sup>) I

## Product Information

### Cat.No.

Kit-1881

### Product Overview

TNKS1 catalyzes the NAD-dependent addition of poly(ADP-ribose) to the substrate proteins.

### Size

384 reactions

### Description

The TNKS1 Histone Ribosylation Assay Kit (Biotin-labeled NAD<sup>+</sup>) is designed to measure Tankyrase 1 (TNKS1) activity for screening and profiling applications. The key to the TNKS1 Histone Ribosylation Assay (Biotin-labeled NAD<sup>+</sup>) is the biotinylated substrate. With this kit, only three simple steps are required for TNKS1 reactions. First, histone proteins are coated on a 384-well plate. Next, the biotinylated TNKS1 substrate is incubated with an assay buffer that contains the TNKS1 enzyme. Finally, the plate is treated with streptavidin-HRP followed by addition of the HRP substrate to produce chemiluminescence that can be measured using a chemiluminescence reader.

### Applications

Great for studying enzyme kinetics and screening small molecular inhibitors for drug discovery and HTS applications.

### Storage

12 months from date of receipt, when stored as directed. Kit components require different storage conditions. Be sure to store each component at the proper temperature upon arrival.

### Kit Components

TNKS1: 2 x 5 µg; -80°C 5 x histone mixture: 2 x 1 ml; -80°C 10x Assay Mixture Containing Biotinylated Substrate: 2 x 300 µl; -80°C 10x PARP Assay Buffer: 2 x 1 ml; -80°C Blocking buffer: 50 ml; +4°C Streptavidin-HRP: 2 x 100 µl; +4°C HRP chemiluminescent substrate A (translucent bottle): 2 x 6 ml; +4°C HRP chemiluminescent substrate B (brown bottle): 2 x 6 ml; +4°C 384-well plate: 1; Room



CREATIVE **BIOMART**<sup>®</sup>  
Assay Kit

## **TNKS1 Histone Ribosylation Assay Kit (Biotin-labeled NAD<sup>+</sup>) I**

Temperature

---

Tel: 1-631-559-9269 1-516-512-3133

Fax: 1-631-938-8127

Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)

45-1 Ramsey Road, Shirley, NY 11967, USA