



## Red Active Caspase-3 Staining Kit

### Product Information

#### Cat

Kit-2372

#### Cat.No.

Kit-2372

### Product Overview

Activation of caspases plays a central role in apoptosis. The Red Active Caspase-3 Staining Kit provides a convenient means for detecting activated caspase-3 in living cells. The assay utilizes the caspase-3 inhibitor DEVD-FMK conjugated to sulfo-rhodamine (Red-DEVD-FMK) as the fluorescent in situ marker. Red-DEVD-FMK is cell permeable, nontoxic, and irreversibly binds to activated caspase-3 in apoptotic cells.

### Applications

Detection of activated caspase-3 in living cells.

### Storage

-20°C

### Shipping

Gel Pack

### Kit Components

Red-VAD-FMK; Wash Buffer; Z-VAD-FMK

### Target Species

Mammalian

**Detection method** Fluorescence microscopy, fluorescence plate reader(Ex/Em = 540/570 nm) and Flow cytometry

### Features & Benefits

Simple one-step procedure; takes only 1-2 hours; Fast and convenient; The red fluorescence label



CREATIVE **BIOMART**<sup>®</sup>  
Assay Kit

## Red Active Caspase-3 Staining Kit

allows for direct detection of activated caspases in apoptotic cells

---

Tel: 1-631-559-9269 1-516-512-3133

Fax: 1-631-938-8127

Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)

45-1 Ramsey Road, Shirley, NY 11967, USA