



Red Active Caspase Staining Kit

Product Information

Cat

Kit-2395

Cat.No.

Kit-2395

Product Overview

Activation of caspases plays a central role in apoptosis. The Red Active Caspase Staining Kit provides a convenient means for detecting activated caspases in; Living cells. The assay utilizes a caspase family inhibitor VAD-FMK conjugated to sulfo-rhodamine (Red-VAD-FMK) as the fluorescent marker. Red-VAD-FMK is cell permeable, nontoxic, and irreversibly binds to activated caspases in apoptotic cells.

Applications

Detection of activated caspases in living cells.

Storage

-20°C

Shipping

Gel Pack

Kit Components

Red-VAD-FMK; Wash Buffer; Z-VAD-FMK

Target Species

Mammalian

Detection method Fluorescence microscopy, fluorescence plate reader(Ex/Em = 540/570 nm) and Flow cytometry

Features & Benefits

Simple one-step procedure; takes only 1-2 hours; Fast and convenient; The red fluorescence label



CREATIVE **BIOMART**[®]
Assay Kit

Red Active Caspase Staining Kit

allows for direct detection of activated caspases in apoptotic cells

Tel: 1-631-559-9269 1-516-512-3133

Fax: 1-631-938-8127

Email: info@creative-biomart.com

45-1 Ramsey Road, Shirley, NY 11967, USA