

Active Native HSV-2 Protein

Cat. No. HSV-2ag-268V Lot. No. (See product label)

SPECIFICATION

| | |
|-------------------------|---|
| Product Overview | Native HSV-2 Protein was produced in BSC-1 cell line. |
| Species | HSV2 |
| Source | BSC-1 |
| Form | Tissue culture medium. |
| Bio-activity | This antigen was inactivated using gamma irradiation. Effectiveness is defined as the inability of an inactivated virus preparation to produce CPE after five days of incubation on a BSC-1 monolayer. |
| Applications | This antigen preparation should be sonicated immediately prior to use to ensure that the preparation is uniform. This preparation may be used as is in a variety of immunoassay formats or may be further purified to meet the needs of a particular assay format. Grade 2 antigen is widely used for both IgG and IgM detection in assays which include EIA with polystyrene and latex solid phases. |
| Stability | The product will be expired one year after receiving if stored properly. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. |
| Warning | This material has been inactivated and should contain no infectious material. However, this material should still be handled with caution as potentially hazardous. |

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA



GENE INFORMATION

Official Symbol HSV-2ag

Synonyms HSV-2ag;

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA