

Acid Phosphatase Activity Colorimetric Assay Kit

Cat. No. Kit-0038 Lot. No. (See product label)

SPECIFICATION

Product Overview	Phosphate concentrations between 1 μ M and 1 mM, with a lower limit of detection of approximately 0.1 nmol, can be directly determined. (500 assays) The assay detects acid phosphatase (AP) activity conveniently in various biological samples.
Size	500 assays
Description	Acid phosphatases (AP) dephosphorylate phosphate groups from phosphate esters under acid conditions. Different acid phosphatase isozymes are found in different organs, and their serum levels are used as a diagnostic for disease in the corresponding organs. Elevated prostatic acid phosphatase levels may indicate the presence of prostate cancer and elevated tartrate-resistant acid phosphatase levels may indicate the bone disease. Acid Phosphatase Assay Kit is a high sensitivity, simple, direct and HTS-ready colorimetric assay designed to measure AP activity in serum and other samples. It is suitable for research and drug discovery. The Acid Phosphatase Activity kit uses p-nitrophenyl phosphate (pNPP) as a phosphatase substrate which turns yellow ($\lambda_{max} = 405$ nm) when dephosphorylated by AP. The kit can detect μ acid phosphatase activity in samples.
Applications	Phosphate concentrations between 1 μ M and 1 mM, with a lower limit of detection of approximately 0.1 nmol, can be directly determined.
Target Species	Mammalian
Usage	This product is furnished for LABORATORY RESEARCH USE ONLY. Not for

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA



diagnostic or therapeutic use.

Kit Components

AP Assay BufferpNPP Substrate (10 TAB)AP EnzymeStop Solution

Detection method


Colorimetric

Compatible Sample Types

Cell Culture SupernatantsPlasmaSerumTissue LysatesUrineCell Lysates

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA