

Caspase-8 Inhibitor Drug Screening Kit (Fluorometric)

Cat. No. Kit-2329 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	<p>Caspases have been shown to play a crucial role in apoptosis induced by various deleterious and physiologic stimuli. Inhibition of caspases can delay apoptosis, implicating a potential role in drug screening efforts. The Caspase-8 Inhibitor Drug Screening Kit provides an effective means for screening caspase inhibitors using fluorometric methods. The assay utilizes synthetic peptide substrate IETD-AFC (AFC, 7-amino-4-trifluoromethyl coumarin). Active caspase-8 cleaves the synthetic substrate to release free AFC which can then be quantified by fluorometry. Compounds to be screened can directly be added to the reaction and the level of inhibition of caspase-8 activity can be determined by comparison of the fluorescence intensity in samples with and without the testing inhibitors.</p>
Applications	Effective means for screening caspase inhibitors using fluorometric methods.
Storage	-20°C
Shipping	Gel Pack
Size	100 assays
Kit Components	2X Reaction Buffer; Caspase Substrate IETD-AFC (1 mM); DTT (1 M); Active Caspase-8 (Lyophilized); Caspase Inhibitor, Z-VAD-FMK (2 mM)
Target Species	Mammalian
Detection method	Fluorescence (Ex/Em 400/505 nm)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA



Features & Benefits

Simple one-step procedure; takes only 1-1.5 hours;

Fast and convenient;

Inhibition of caspase-8 activity can be determined by comparison of the fluorescence intensity in samples with and without the testing inhibitors.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA