

## Glucose-6-Phosphate Assay Kit

Kit-2490

Lot. No. (See product label)

### Specification

#### Product Overview

Glucose-6-Phosphate (g6p) is glucose sugar phosphorylated on carbon 6. Most of the glucose entering cells is phosphorylated to G6P. G6P has three primary fates within the cell. It lies at the start of two major metabolic pathways: glycolysis and the pentose phosphate pathway. In addition to these metabolic pathways, glucose 6-phosphate may also be converted to glycogen or starch for storage. The Glucose-6-Phosphate Assay Kit provides a simple, and automation-ready procedure for measuring G6P concentration. G6P is oxidized by glucose-6-phosphate dehydrogenase and the formed NADPH is coupled to the formazan (WST-8) chromogen. The intensity of the product color, measured at 460 nm, is proportional to the G6P concentration in the sample.

#### Storage

-20°C

#### Shipping

On Ice

#### Size

100 tests

#### Detection method

OD460nm

#### Compatible Sample Types

Plasma, serum, tissue and culture media etc

#### Features & Benefits

Fast and sensitive. Linear detection range: 10 to 1000  $\mu$ M G6P.  
Convenient and high-throughput. Homogeneous 'mix-incubate-measure' type assay. Can be readily automated on HTS liquid handling systems for processing thousands of samples per day.

#### Assay time

20 min

#### Sensitivity

10  $\mu$ M

For Research Use Only

Creative BioMart. All rights reserved

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: +1-631-559-9269 Fax: +1-631-938-8127

E-mail: [info@creative-biomart.com](mailto:info@creative-biomart.com)

[www.creativebiomart.net](http://www.creativebiomart.net)