

Protein A Beads

Cat. No. Protein A-28A **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Protein A beads may be used for purification of IgG subclass of antibodies from several species. For a complete use guide, please follow the protocol section. When used and stored properly, the beads may be used repeatedly for up to 100 times.
Cat.no	Protein A-28A
Matrix	Highly cross-linked agarose, 4%
Functional group	Protein A
Average particle size	45-165 µm
Ligand concentration	1-2 mg/ml
Dynamic binding capacity	4-7 mg IgG per ml
Recommended flow rate	1-3 ml/min/ cm ²
pH working range	3-11
Elution buffer	2% sucrose, adjusted pH to 2-3 by CH ₃ COOH.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA