

## Native Tyrophagus Longior Extract Reference Standard

**Cat. No.** Tyrophagus Longior Extract-33ARS    **Lot. No.** (See product label)

### SPECIFICATION

<b>Species</b>	Tyrophagus longior mite
<b>Source</b>	Tyrophagus longior mite
<b>Tag</b>	Non
<b>Form</b>	Freeze-dried powder
<b>Stability</b>	Accelerated degradation studies have shown freeze dried allergen extracts stored in unopened ampoules at -20°C to be extremely stable over a number of years. Reference Materials should be stored on receipt as indicated on the label.
<b>Storage</b>	Unopened ampoules should be stored at -20°C.
<b>contents</b>	The extract was thawed at room temperature, filtered through a sterile HA Millipore, stored overnight at 4°C and distributed in 1.0ml volumes at room temperature into approximately 400 ampoules.
<b>Shipping</b>	Because of the inherent stability of lyophilized material, may ship these materials at ambient temperature.
<b>USE OF MATERIAL</b>	No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution. The total contents of the ampoule should be reconstituted with 1.0ml distilled water and dissolved gently by swirling to avoid froth. No attempt should be made to weigh

 Tel: 1-631-559-9269    1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)     Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA



out any portion of the freeze dried material. The reconstituted reagent should be used as soon as possible after reconstitution. The standard can also be used in in vitro assays such as the radioallergosorbent test (RAST), RAST inhibition, quantitative immunoelectrophoresis (CIE & CRIE), SDS-PAGE and immunoblotting.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA