

Direct Glucose Uptake Assay Kit

Cat. No. Kit-2229 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Glucose uptake is one of the key processes for cellular glucose metabolism. The study of glucose uptake can provide important information for understanding glucose metabolism and regulation in normal and disease development such as diabetes. EZCell™ Glucose Uptake Assay kit is simple, ultra-sensitive and easy to use. A specific hexokinase inhibitor that inhibits hexokinase, the first enzyme metabolizing glucose in cells is used to arrest glucose consumption after its uptake. Glucose Uptake is measured by using a set of enzymatic reactions that specifically oxidize glucose producing intermediates that react with the OxiRed; Probe generating a fluorescence signal (Ex/Em=535/587 nm). The fluorescence signal is directly proportional to the amount of glucose that has been taken up and accumulated inside the cells. Unlike other kits detecting glucose derivatives, this glucose uptake assay provides a direct, powerful tool for studying this process as well as for screening and characterization of drugs that regulate glucose uptake during normal and disease development.

Applications

Simple procedure; Fast and convenient; The assay is easy, non-radioactive, and highly sensitive.

Storage

-20°C

Shipping

Gel Pack

Size

50 assays

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Kit Components	Assay Buffer; OxiRed™ Probe (in DMSO); Enzyme Mix; Hexokinase Inhibitor; Glucose Standard (100 mM); Glucose (1 M, Sterile)
Target Species	Mammalian
Detection method	Fluorescence (535/587 nm)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA