

Glucose Uptake Fluorometric Assay Kit

Cat. No. Kit-2250 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Glucose uptake is an important biological process for studying cell signaling and glucose metabolism. Among many different methods available for measuring glucose uptake, 2-deoxyglucose (2-DG) has been widely used because of its structural similarity to glucose. As with glucose, 2-DG can be taken up by glucose transporters, and metabolized to 2-DG-6-phosphate (2-DG6P). 2-DG6P, however, cannot be further metabolized, and thus accumulates in the cells. The accumulated 2-DG6P is, therefore, directly proportional to 2-DG (or glucose) uptake by cells. In glucose uptake fluorometric assay kit, the accumulated 2-DG6P is enzymatically oxidized and coupled to PicoProbe, which generates fluorescence in the presence of NADPH. This easy to use non-radioactive kit can detect glucose uptake as low as 50 pmol/well and can be used for a variety of cell types.

Applications

This easy to use non-radioactive kit can detect glucose uptake as low as 50 pmol/well and can be used for a variety of cell types.

Storage

-20°C

Shipping

Gel Pack

Size

100 assays

Kit Components

Extraction Buffer; Neutralization Buffer; 2-Deoxyglucose (2-DG, 10 mM); 2-DG Uptake Assay Buffer; Enzyme Mix (Lyophilized); (in DMSO); 2-DG6P Standard (Lyophilized)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Target Species	Mammalian
Detection method	Fluorescence(Ex/Em = 535/587 nm)
Features & Benefits	Simple procedure; Fast and convenient; The assay is easy, non-radioactive, and highly sensitive

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA