

Kit for Saliva Exosomes

Cat. No. Kit-2308 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Exosome isolation from biofluids or cell culture media and FACS analysis of exosome markers.
Applications	Kits ensure exosome isolation and exosome marker characterization via FACS, Exosome comprehensive profiling, User friendly and suitable for multiple marker analyses, No initial exosome purification required, Lyophilized Exosome Standards for positive control included.
Storage	4°C
Shipping	Gel Pack
Size	20 reactions
Kit Components	ExoPure Reagent; FACS-beads; Primary Antibody; Secondary Antibody; Sample Buffer (1X); Exosome Standards (100 µg)
Target Species	Saliva
Detection method	FACS Analysis
Features & Benefits	Exosome isolation and exosome marker characterization via FACS

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA