

## Monoamine Oxidase B (MAO-B) Inhibitor Screening Kit (Fluorometric)

Cat. No. Kit-2362    Lot. No. (See product label)

### SPECIFICATION

#### Product Overview

Monoamine oxidases (MAO, EC 1.4.3.4) are a family of enzymes that can oxidize a wide variety of endogenous primary amines. Two isoforms, MAO-A and MAO-B, have been identified based on their substrate, inhibitor specificity, and tissue localization. MAO-B can oxidize primary amines, but its list of specific substrates (i.e. benzylamine, phenylethylamine) is more limited compared to MAO-A. MAO-B is a mitochondrial-bound enzyme that is ubiquitously expressed throughout the brain and other tissues. It has been investigated in numerous studies including Parkinson's disease, Alzheimer's, and tobacco addiction. Specific MAO-B inhibitors such as selegiline, & rasagiline have been used to treat Parkinson's patients, but their benefits are considered rather modest. MAO-B Inhibitor Screening Kit offers a rapid, simple, sensitive, and reliable test suitable for high-throughput screening of MAO-B inhibitors. The assay is based on the fluorometric detection of H<sub>2</sub>O<sub>2</sub>, one of the byproducts generated during the oxidative deamination of MAO substrate.

#### Applications

Screening/studying/characterizing MAO-B inhibitors

#### Storage

-20°C

#### Shipping

Gel Pack

#### Size

100 assays

#### Kit Components

MAO-B Assay Buffer; OxiRed; Probe (in DMSO); MAO-B Enzyme (Lyophilized);

 Tel: 1-631-559-9269    1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)     Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

|                                |                                                                                                                                                            |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                | MAO-B Substrate (Lyophilized); Developer (Lyophilized); Inhibitor Control (Clorgyline) (Lyophilized)                                                       |
| <b>Detection method</b>        | Fluorescence (Ex/Em = 535/587 nm)                                                                                                                          |
| <b>Features &amp; Benefits</b> | Rapid & sensitive; High-throughput assay to screen/study/characterize potential inhibitors of monoamine oxidase B; Includes Inhibitor Control [Selegiline] |

 Tel: 1-631-559-9269 1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA