

Red Active Caspase-3 Staining Kit

Cat. No. Kit-2372 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Activation of caspases plays a central role in apoptosis. The Red Active Caspase-3 Staining Kit provides a convenient means for detecting activated caspase-3 in; Living cells. The assay utilizes the caspase-3 inhibitor DEVD-FMK conjugated to sulfo-rhodamine (Red-DEVD-FMK) as the fluorescent in situ marker. Red-DEVD-FMK is cell permeable, nontoxic, and irreversibly binds to activated caspase-3 in apoptotic cells.
Applications	Detection of activated caspase-3 in living cells.
Storage	-20°C
Shipping	Gel Pack
Kit Components	Red-VAD-FMK; Wash Buffer; Z-VAD-FMK
Target Species	Mammalian
Detection method	Fluorescence microscopy, fluorescence plate reader(Ex/Em = 540/570 nm) and Flow cytometry
Features & Benefits	Simple one-step procedure; takes only 1-2 hours; Fast and convenient; The red fluorescence label allows for direct detection of activated caspases in apoptotic cells

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA