

Red Active Caspase-8 Staining Kit

Cat. No. Kit-2373 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Activation of caspases plays a central role in apoptosis. The Red Active Caspase-8 Staining Kit provides a convenient means for detecting activated caspase-8 in; Living cells. The assay utilizes a caspase-8 inhibitor IETD-FMK conjugated to sulfo-rhodamine (Red-IETD-FMK) as the fluorescent marker. Red-IETD-FMK is cell permeable, nontoxic, and irreversibly binds to activated caspase-8 in apoptotic cells.
Applications	Detection of activated caspase-8 in living cells.
Storage	-20°C
Shipping	Gel Pack
Kit Components	Red-IETD-FMK; Wash Buffer; Z-VAD-FMK
Target Species	Mammalian
Detection method	Fluorescence microscopy, fluorescence plate reader(Ex/Em = 540/570 nm) and Flow cytometry
Features & Benefits	Simple one-step procedure; takes only 1-2 hours; Fast and convenient; The fluorescence label allows detection of activated caspase-8 in apoptotic cells

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA