

Annexin V-FITC Apoptosis Kit

Cat. No. Kit-2376 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Annexin V Apoptosis Detection Kit is based on the observation that soon after initiating apoptosis, cells translocate the membrane phosphatidylserine (PS) from the inner face of the plasma membrane to the cell surface. Once on the cell surface, PS can be easily detected by staining with a fluorescent conjugate of Annexin V, a protein that has a high affinity for PS. The one-step staining procedure takes only 10 minutes. Detection can be analyzed by flow cytometry or by fluorescence microscopy. The kit can differentiate apoptosis vs necrosis when performing both Annexin V-FITC and PI staining.
Applications	Detect early/middle stages of apoptosis; differentiate apoptosis from necrosis.
Storage	4°C
Shipping	Gel Pack
Kit Components	Annexin V-FITC; 1X Binding Buffer; Propidium Iodide (PI)
Target Species	Mammalian
Detection method	Flow cytometry(Ex = 488 nm; Em = 530 nm) and fluorescence microscopy
Features & Benefits	Simple one step staining procedure in 10 minutes; Fast and convenient; Kit can differentiate apoptosis vs necrosis when performing both Annexin V-FITC and PI staining

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA