

Fluorescein Active Caspase-8 Staining Kit

Cat. No. Kit-2384 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Activation of caspases plays a central role in apoptosis. The CaspGLOW;Fluorescein Active Caspase-8 Staining Kit provides a convenient means for sensitive detection of activated caspase-8 in; Living cells. The assay utilizes the caspase-8 inhibitor, IETD-FMK, conjugated to FITC (FITC-IETD-FMK) as a fluorescent in situ marker. FITC-IETD-FMK is cell permeable, nontoxic, and irreversibly binds to activated caspase-8 in apoptotic cells.
Applications	Sensitive detection of activated caspase-8 in living cells.
Storage	-20°C
Shipping	Gel Pack
Kit Components	FITC-ATAD-FMK; Wash Buffer; Z-VAD-FMK
Target Species	Mammalian
Detection method	Fluorescence microscopy, Flow cytometry and fluorescence plate reader (Ex. = 485 nm and Em. = 535 nm)
Features & Benefits	Simple one-step procedure; takes only 1-2 hours; Fast and convenient; The FITC label allows for direct detection of the activated caspases in apoptotic cells

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA