

NADP/NADPH Assay Kit

Cat. No. Kit-2434 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Pyridine nucleotides play an important role in metabolism and, thus, there is continual interest in monitoring their concentration levels. Quantitative determination of NADP⁺/NADPH has applications in research pertaining to energy transformation and redox state of cells or tissue. Simple, direct and automation-ready procedures for measuring NADP⁺/NADPH concentration are very desirable. NADP⁺/NADPH assay kit is based on a glucose dehydrogenase cycling reaction, in which the formed NADPH reduces a formazan (MTT) reagent. The intensity of the reduced product color, measured at 565 nm, is proportionate to the NADP⁺/NADPH concentration in the sample. This assay is highly specific for NADP⁺/NADPH and is not interfered by NAD⁺/NADH. Our assay is a convenient method to measure NADP, NADPH and their ratio.

Storage

-20°C

Shipping

On Ice

Size

100 tests

Detection method

OD565nm

Compatible Sample Types

Cell, tissue extracts etc

Features & Benefits

Sensitive and accurate. Detection limit 0.1 µM, linearity up to 10 µM NADP⁺/NADPH in 96-well plate assay.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Convenient. The procedure involves adding a single working reagent, and reading the optical density at time zero and 30 min at room temperature. No 37°C heater is required.

High-throughput. Can be readily automated as a high-throughput 96-well plate assay for thousands of samples per day.

Assay time	30 min
Sensitivity	0.1 μ M

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA