

D-Mannitol Assay Kit

Cat. No. Kit-2460 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	D-MANNITOL (C ₆ H ₁₄ O ₆) is a sugar alcohol used in dietary supplement, sweetener, intestinal permeability test for leaky gut, etc. It also serves as a coating for hard candies, dried fruits, and chewing gums due to its low ability to attract and hold water molecules. In addition, it is an osmoprotectant for plants and is used clinically in osmotherapy to reduce intracranial pressure. D-mannitol assay kit is based on mannitol dehydrogenase catalyzed oxidation of D-mannitol, which generates D-fructose and NADH that reduces a formazan (MTT) dye. The intensity of product color, measured at 565 nm is directly proportional to D-mannitol concentration in the sample.
Storage	-20°C
Shipping	On Ice
Size	100 tests
Detection method	OD565nm
Compatible Sample Types	Biological, food, beverage and agricultural products
Features & Benefits	Fast and sensitive. Use of 20 µL sample. Linear detection range 0.007 to 3 mM mannitol in 96-well plate assay. Convenient. The procedure involves adding a single working reagent, and reading the absorbance after 30 minutes. Room temperature assay. No 37°C heater is needed.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

High-throughput. "Add-mix-read" type assay. Can be readily automated as a high-throughput 96-well plate assay for thousands of samples per day.

Assay time	30 min
-------------------	--------

Sensitivity	7 μ M
--------------------	-----------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA