

Glucose-6-Phosphate Assay Kit

Cat. No. Kit-2490 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Glucose-6-Phosphate (g6p) is glucose sugar phosphorylated on carbon 6. Most of the glucose entering cells is phosphorylated to G6P. G6P has three primary fates within the cell. It lies at the start of two major metabolic pathways: glycolysis and the pentose phosphate pathway. In addition to these metabolic pathways, glucose 6-phosphate may also be converted to glycogen or starch for storage. The Glucose-6-Phosphate Assay Kit provides a simple, and automation-ready procedure for measuring G6P concentration. G6P is oxidized by glucose-6-phosphate dehydrogenase and the formed NADPH is coupled to the formazan (WST-8) chromogen. The intensity of the product color, measured at 460 nm, is proportional to the G6P concentration in the sample.
-------------------------	---

Storage	-20°C
----------------	-------

Shipping	On Ice
-----------------	--------

Size	100 tests
-------------	-----------

Detection method	OD460nm
-------------------------	---------

Compatible Sample Types	Plasma, serum, tissue and culture media etc
--------------------------------	---

Features & Benefits	Fast and sensitive. Linear detection range: 10 to 1000 µM G6P. Convenient and high-throughput. Homogeneous "mix-incubate-measure" type assay. Can be readily automated on HTS liquid handling systems for processing thousands
--------------------------------	--

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

of samples per day.

Assay time	20 min
-------------------	--------

Sensitivity	10 uM
--------------------	-------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA