

Recombinant Full Length Vibrio Cholerae Serotype O1 Type 4 Prepilin-Like Proteins Leader Peptide-Processing Enzyme(Tcpj) Protein, His-Tagged

Cat. No. RFL10841VF **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant Full Length Vibrio cholerae serotype O1 Type 4 prepilin-like proteins leader peptide-processing enzyme(tcpJ) Protein (A5F385) (1-253aa), fused to N-terminal His tag, was expressed in E. coli.
Species	Vibrio cholerae serotype O1
Source	E.coli
ProteinLength	Full Length (1-253)
Form	Lyophilized powder
AA Sequence	MEYVYLILFSIVSLILGSFSNVVIYRLPRKILLKNHFFYDIDSNRSMCPKCGNKISWYDN VPLLSYLLLHGKCRHCDEKISLSYFIVELSFIIAFPIYWLSTDWVDSFVLLGLYFILFN L FVIDFKSMLLPNLLTYPIFMLAFIYVQNPALTVESSIIGGFAAFIISYVSNFIVRLFK RID VMGGGDIKLYTAIGTLIGVEFVPYLFLLSSIIAFIHWFFARVSCRYCLYIPLGPSII ISFVI VFFSIRLM
Purity	Greater than 90% as determined by SDS-PAGE.
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.
Storage Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
GENE INFORMATION	
Gene Name	tcpJ
Synonyms	tcpJ; VC0395_A0364; VC395_0855; Prepilin leader peptidase/N-methyltransferase [Includes: Leader peptidase; Prepilin peptidase; N-methyltransferase;]
UniProt ID	A5F385

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA