

Recombinant Full Length Botryotinia Fuckeliana Golgi Apparatus Membrane Protein Tvp38(Tvp38) Protein, His-Tagged

Cat. No. RFL10681BF **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant Full Length Botryotinia fuckeliana Golgi apparatus membrane protein tvp38(tvp38) Protein (A6RIB9) (1-392aa), fused to N-terminal His tag, was expressed in E. coli.
Species	Botryotinia fuckeliana
Source	E.coli
ProteinLength	Full Length (1-392)
Form	Lyophilized powder
AA Sequence	MPADYSSTARALALPISPLETPSSPDQHVERPPWSRRISSNVRRANNSASPYSNQP TTSF KDQILRKAELQQRKFVSTWQKLSPLQKCLAVGAAVLNIIILALFLVYQHQIFASL APFAE RWRDMRGGWMILWAMTFVAAFPLIGYSSTITIAGFVYGVPKGWAIVASAT VAGSLCSFL ASRTILSSYVHRLVGKDKRFEALATLKHDGIKILCMIRLCPLPYLSNA AVATFPTVHP LNYALATALVTPKLFHVFIGSRLGSLAGDEEMDASTKLINYASIIIGAG LGATVGYVIY QRTMARAKELEIEELEANGDVAAGRVA AEYSDANNDDAALMNDD DISLWDNEDAQTGY TDFVDEDEADV FASGDLDEEGNIGGKGKTAGT
Purity	Greater than 90% as determined by SDS-PAGE.
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

	for up to one week.
Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.
Storage Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

GENE INFORMATION

Gene Name	tpv38
Synonyms	tpv38; BC1G_00190; Golgi apparatus membrane protein tpv38
UniProt ID	A6RIB9

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA