

Active Recombinant Heat-labile UDG

Cat. No. UDG-42 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Heat-labile UDG cloned from psychrophilic marine bacterium and expressed in E. coli.
Species	Psychrophilic Marine Bacterium
Source	E.coli
Description	Heat-labile UDG (Uracil DNA Glycosylase) can catalyze the hydrolysis of the N-glycosidic bond from deoxy uridine to release uracil in ssDNA and dsDNA. This product is sensitive to heat and the enzyme activity is irreversibly inactivated at temperature above 50 centigrade. It is suitable for PCR, qPCR, RT-PCR and RT-qPCR. It can also be used to remove aerosol contamination from dU-containing PCR products.
Bio-activity	The specific activity of Heat-labile UDG was determined to be 1 Unit/ μ L.
Unit Definition	One unit (U) is defined as the amount of enzyme that releases 1 nmol of uracil from dU-containing DNA within 1h at 37 centigrade in the reaction system containing 70 mM of Tris-HCl, pH 7.5, 10 mM of NaCl, 1 mM of EDTA, 100 μ g/ml of BSA.
Stability	One year at -20 centigrade from date of shipment.
Storage	Store product at -20 centigrade. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For best performance, avoid repeated handling and multiple freeze/thaw cycles.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Storage Buffer	Recombinant protein stored in 20 mM Tris-HCl, pH 8.0, 0.1 mM EDTA, 100 mM KCl, 1 mM DTT, 50% glycerol, 0.5% NP-40, 0.5% Tween-20.
-----------------------	---

Shipping	Ice pack
-----------------	----------

GENE INFORMATION

Official Symbol	UDG
------------------------	-----

Synonyms	Uracil DNA Glycosylase; Heat-labile UDG; UDG
-----------------	--

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA