

CD45R protein-coupled magnetic MicroBeads

Cat. No. CD45R-238M **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Species Rat

Capacity for 1×10^9 total cells


Background

The CD45R antigen is expressed throughout B cell development from the pro-B cell stage on, but it is down-regulated during the differentiation of B cells into plasma cells. Expression levels of CD45R are different among distinct B cell subpopulations. Memory B cells are CD45R^{low/-}. CD45R is also weakly expressed on a small subset of peripheral T cells, but is not detected on thymocytes.

Application

B cells isolated from bone marrow or peripheral lymphoid organs can be used for in vitro or ex vivo analysis of B cell lymphopoiesis and B cell migration. B cells purified with CD45R MicroBeads can be used for studies on the mechanism and regulation of immunoglobulin gene rearrangement.

Analysis



 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Isolation of B cells from rat spleen using CD45R MicroBeads, a MiniMACS Separator, and an MS Column.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA