

Dextrin Sepharose High Performance MBP-tagged protein purification resin

Cat. No. Dextrin-01HP **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Dextrin Sepharose High Performance is a robust and stable affinity resin for purification of proteins tagged with maltose binding protein (MBP). Fast, one-step purification of MBP-tagged proteins. High-resolution purification and elution in narrow peaks, minimizing the need for further concentration steps. Can be used for repeated purification: can be easily regenerated using 0.5 M NaOH. Physiological conditions and mild elution preserve target protein activity.
Chromatography technique	Affinity Chromatography
Chemical stability	Stable in commonly used aqueous buffers, 0.5 M NaOH (regeneration and cleaning)
pH stability, CIP	2–13
Storage	2 to 8°C, 20% Ethanol
Binding	Approx. 7 mg MBP- Δ Sal/ml resin Approx. 16 mg MBP- β Gal/ml resin
Size	25ml; 100ml
Matrix	Highly cross-linked agarose, 6%
Average particle size	~34 μ m

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Ligand	Dextrin
Pack size	25 mL
Flow velocity	< < 300 cm/h

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA