

Buffy Coat CD4 protein-coupled magnetic MicroBeads

Cat. No. CD4-4832M **Lot. No.** (See product label)

SPECIFICATION

Product Overview


Species Human

Capacity for 1 buffy coat

Background The CD4 antigen is expressed on T helper cells and at a lower level on monocytes and dendritic cells. The CD4 antigen is an accessory molecule in the recognition of MHC class II/peptide complexes by the TCR heterodimer on CD4+ T helper cells.

Application The purified CD4+ cells are well suited for further flow cytometric, functional, or molecular analysis.

Analysis



Separation of a buffy coat sample using the StraightFrom Buffy Coat CD4 MicroBead Kit and the MultiMACS Cell24 Separator Plus with the Single-Column Adapter and Whole Blood Columns. Cells were fluorescently stained with CD4-PE, CD14-APC, as well as CD45-VioBlue and analyzed by flow cytometry using the MACSQuant Analyzer. Cells were triggered via CD45-VioBlue, cell debris and dead cells were excluded from the analysis based on scatter signals and propidium iodide

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

fluorescence.

Tel: 1-631-559-9269 1-516-512-3133

Email: info@creative-biomart.com Fax: 1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA