

Leukopak REAlease CD56 protein-coupled magnetic MicroBeads

Cat. No. CD56-6832M **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Species Human

Capacity for 0.5 Leukopaks

Background

The StraightFrom Leukopak REAlease MicroBead Kit was developed for the rapid positive isolation of CD56+ cells directly from Leukopak using the MultiMACS Cell 24 Separator Plus for semi-automated separation. No sample preparation is required, including density gradient centrifugation, erythrocyte lysis, or cell count. The Multi-24 Column Block is included in the kit to allow direct processing of a 1/2 Leukopak.

Application

CD56+ cells isolated using StraightFrom Leukopak REAlease CD56 MicroBeads are ideally suited for further magnetic separation, i.e. for the depletion of CD56+CD3+ cells with CD3 MicroBeads. The cells are suited for further flow cytometric analyses, molecular biology applications and functional studies.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA