

Spleen CD19 protein-coupled magnetic MicroBeads

Cat. No. CD19-2913M **Lot. No.** (See product label)

SPECIFICATION


Product Overview

Species Mouse

Capacity for 40 spleens

Background

The StraightFrom Spleen CD19 MicroBead Kit, mouse has been developed for the positive selection of CD19⁺ cells directly from mouse spleen by using the gentleMACS Dissociators in combination with enzymatic digestion and MACS Technology.



Analysis

Separation of a spleen sample using the StraightFrom Spleen CD19 MicroBead Kit, mouse and the autoMACS NEO Separator. Cells were fluorescently stained with CD19-PE and CD45R (B220)-APC-Vio 770 and analyzed by flow cytometry using the MACSQuant Analyzer. Cell debris and dead cells were excluded from the analysis based on scatter signals and 7-AAD fluorescence.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA