

Active TACSTD2 protein-coupled magnetic Beads

Cat. No. TACSTD2-4932M **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Species	Human
Source	HEK293
ProteinLength	Gln 31 - Thr 274
Appearance	Powder mixture
Particle size	2 µm
Bio-activity	Immobilized 32.6 µg TROP-2/1 mg beads can bind Anti-TROP-2 mAb with an EC50 of 0.8505 µg/mL.
Coupled amount of protein	≈1032 pmol (32.6 µg) TROP-2/mg beads
Capacity	>266 pmol (40µg) Anti-trop-2 antibody/mg beads
Formulation	PBS, pH7.4, with 10% Trehalose
Reconstitution	2 mL sterile deionized water (1 mg beads/mL)
Background	The Human TROP-2-coupled Magnetic Beads is produced by coupling biotinylated TROP-2 heterodimer to streptavidin-conjugated magnetic beads. Because

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Streptavidin (SA) has an extraordinarily high affinity for biotin with a dissociation constant (K_d) on the order of 10^{-14} mol/L, the biotinylated protein can bind to the SA beads irreversibly. The pre-coupled beads are ready to use for capturing anti-TROP-2 antibody from your sample with high specificity.

Application

This product is intended for immunocapture, biopanning and flow cytometry. This product is produced non-sterile.

Storage

Upon receipt, please store the Beads at -20°C for 1 year in lyophilized state.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA