

Recombinant allergen Ole e 1 for Olea europaea (Olive tree)

Cat. No. IVP-016 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Ole e 1 is the most prevalent allergen of olive tree pollen (<i>Olea europaea</i>). This allergen has been prepared as a recombinant mature protein fused to a his-tag.
Source	P.pastoris
Molecular Mass	Determined by SDS-PAGE, the protein smear is located between molecular markers of 35,000 and 25,000 Da, while relative molecular mass calculated from amino acid sequence, without glycosylation, is approx. 18,887.36 Da.
Purity	>95% in a SDS-PAGE
Incidence	The tested incidence of the protein is 100%, in front of 70% according to the bibliography.
IVD	Serodiagnosis of allergy to Olive tree
Stability	The protein will remain stable for four years if stored frozen.
Observation	Proteins should be maintained frozen at high concentrations. In order to defrost the protein, maintain the aliquot at 25°C without shaking to avoid aggregation. Prior making test dilutions and after defrosting the protein, is recommended to remove possible protein aggregates by centrifuging the stock solution, avoiding alterations in the immobilization of the biomolecule to the solid surface.
Storage	Protein is shipped with dry ice. Upon arrival, it should be aliquoted in order to avoid

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

repeated freezing and thawing cycles and stored at -20°C to -80°C.

Concentration

1.20mg/ml

Notes

During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 ul or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the containers cap.

Although recombinant antigens are expressed in non-pathogenic P.pastori and bacterial integrity is destroyed during purification, the antigen preparation should be handled as potentially infectious.

Type

Recombinant

ClassID 1

Proteins similar to Ole e 1

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA