

MK2 CoIP Agarose kit

Cat. No. MK2-02A **Lot. No.** (See product label)

SPECIFICATION


Product Overview	The MK2-Trap Agarose Kit contains MK2-Trap Agarose, lysis, wash, and elution buffers for efficient Immunoprecipitation of MK2. • Fast, reliable & efficient one-step Immunoprecipitation • Ready-to-use • No heavy & light antibody chains on the SDS-PAGE • Stable under harsh washing conditions • Short incubation time (1 h at 4°C)
Specificity/Target	MK2 (MAP kinase-activated protein kinase 2, MAPKAP-K2)
Tested Reactivity	Human, Mouse, Hamster Posttranslational modifications: MK2-Trap recognizes unphosphorylated MK2 and Phospho-MK2 (Thr222). Specificity on Phospho-MK2 (Thr334) was not tested.
Conjugate	Agarose beads; bead size: ~ 90 µm (cross-linked 4 % agarose beads)
Elution buffer	SDS sample buffer 0.2 M glycine pH 2.5
Type	Nanobody
Storage Buffer	20% ethanol
Applications	IP, CoIP

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Western Blot



Tel: 1-631-559-9269 1-516-512-3133

Email: info@creative-biomart.com Fax: 1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA