

Native Porcine Gastrin II Reference Standard

Cat. No. Gastrin II-22PRS **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Species	Porcine
Source	Porcine
Tag	Non
Form	Freeze-dried powder
Stability	Reference Materials should be stored on receipt as indicated on the label. In addition, once reconstituted, diluted or aliquoted, users should determine the stability of the material according to their own method of preparation, storage and use.
Usage	Calibrant for porcine gastrin II assays
Storage	Store at -20°C or below in the dark.
contents	Each ampoule contains approximately: Highly-purified porcine gastrin II 10µg Sucrose 5mg The residue of approx. 0.1 M Sorensen's buffer
Shipping	Because of the inherent stability of lyophilized material, may ship these materials at ambient temperature.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

UNITAGE	Each ampoule of the Standard contains 10 units (by definition)
USE OF MATERIAL	No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution. For all practical purposes each ampoule contains the same amount of the same materials. Dissolve the total contents in a known amount of suitable buffer solution with carrier protein (free of peptidase), where extensive dilution is required, to minimise loss by surface adsorption. For economy of use it is recommended that the solution be subdivided into several small containers, which are frozen rapidly to below -40°C and then stored below -40°C in the dark. Repeated freezing and thawing should be avoided.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA