

Hep2 Whole Cell Lysate

Cat. No. Hep2-01HL **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Hep2 Whole Cell Lysate
Species	Human
Source	HEp2 (HeLa contaminant)
Lysate preparation	The cells were grown in MEM supplemented with 10% FBS (Fetal Bovine Serum). The lysate was prepared by first washing the cells in PBS supplemented with EDTA and a cocktail of protease inhibitors (see below). Washed cells are then incubated on ice in modified RIPA buffer containing 150 mM sodium chloride, 50 mM Tris Cl, pH 7.4, 1 mM EDTA, 1.0% NP-40, 0.25% sodium deoxycholic acid to lyse the cells. Protein integrity is ensured using a cocktail of protease inhibitors with broad specificity for the inhibition of aspartic, cysteine, and serine proteases as well as aminopeptidases (0.5 mM AEBSF HCl, 0.4 mM Aprotinin, 25 µM Bestatin, 7.5 µM E-64, 10 µM Leupeptin Hemisulfate and 5 µM Pepstatin A). Cell debris was removed by centrifugation. Protein concentration was determined by biuret assay using a commercially available kit. The protein concentration was adjusted to 2 mg/mL and then an equal volume of 2× SDS-PAGE sample buffer was added.
Recommended Usage	Ready-to-use lysates are especially prepared as positive controls for separation by SDS-PAGE and subsequent western blot analysis. Lysates are prepared in denaturing buffer WITHOUT dissociating agents (i.e. no 2- mercaptoethanol or dithiothreitol has been added).

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Storage Store at 2-8 centigrade for continuous use. For extended storage, freeze working aliquots at -70 centigrade. Repeated freezing and thawing is not recommended. Under proper storage conditions the shelf life is half a year from the date of receipt.

Storage Buffer 1X SDS-PAGE Sample Buffer (62.5 mM Tris HCl, 2% SDS, 10% Glycerol and 0.005% bromophenol blue, pH 6.8)

Concentration 2 mg/mL

SDS-PAGE



 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA