

## Pneumococcal Type 9V Antigen, CRM197 Conjugated

**Cat. No.** Pneumococcal-47B    **Lot. No.** (See product label)

### SPECIFICATION

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b>  | CRM197 Conjugated Pneumococcal Type 9V Antigen                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Species</b>           | Pneumococcal                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Source</b>            | Pneumococcal                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Conjugation/Label</b> | CRM197                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Description</b>       | In most cases, once pneumococcal pneumonia has been identified, doctors will prescribe antibiotics. These antibiotic usually help alleviate and eliminate symptoms between 12 and 36 hours after being taken. Despite most antibiotics' effectiveness in treating the disease, sometimes the bacteria can resist the antibiotics, causing symptoms to worsen. Additionally, age and health of the infected patient can contribute to the effectiveness of the antibiotics. |
| <b>Purification</b>      | IEC & HPLC                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Purity</b>            | > 95%                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Application</b>       | Immunogen, WB, ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Storage</b>           | Store the lyophilized powder at 4 centigrade for 6 months. It is recommended to aliquote to ~10ul per vial and store at -20 centigrade for long term storage. Avoid repeated freezing and thawing cycles.                                                                                                                                                                                                                                                                  |

 Tel: 1-631-559-9269    1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)     Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

|                       |                            |
|-----------------------|----------------------------|
| <b>Storage Buffer</b> | Powder or 0.1M PBS, PH 7.0 |
|-----------------------|----------------------------|

## GENE INFORMATION

|                 |              |
|-----------------|--------------|
| <b>Synonyms</b> | Pneumococcal |
|-----------------|--------------|

 Tel: 1-631-559-9269 1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA