

Recombinant Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720) flag protein, His-tagged

Cat. No. flag-0465S **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720) flag protein, His-tagged(P06179)(1-495 aa), fused with His tag, was expressed in E. coli.
Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Source	E.coli
ProteinLength	1-495 aa
Form	Lyophilized from a 0.2 µm filtered concentrated solution in 1×PBS, 3 % Trehalose, 0.15 % Tween-20.
Molecular Mass	Approximately 52.7 kDa, a single non-glycosylated polypeptide chain containing 503 amino acids, with Leu, Glu and 6 × His at C-terminus.
AASequence	MAQVINTNSL SLLTQNNLNK SQSALGTAIE RLSSGLRINS AKDDAAGQAI ANRFTANIKG LTQASRNAND GISIAQTTEG ALNEINNNLQ RVRELAVQSA NSTNSQSDLD SIQAEITQRL NEIDRVSGQT QFNGVKVLAQ DNTLTIQVGA NDGETIDIDL KQINSQTLGL DTLNVQQKYK VSDTAATVTG YADTTIALDN STFKASATGL GGTDQKIDGD LKFD DTTGKY YAKVTVTGGT GKDGYYEVS DKTNGEVTLA GGATSPLTGG LPATATEDVK NVQVANADLT EAKAALTAAG VTGTASVVKM SYTDNNGKTI DGGLAVKVGD DYYSATQNKD GSISINTTKY TADDGTSKTA LNKLGADGK TEVVSIGGKT YAASKAEGHN FKAQPDLEA AATTENPLQ KIDAALAQVD TLRSDLGAVQ NRFNSAITNL GNTVNNLTSA

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

	RSRIEDSDYA TEVSNMSRAQ ILQQAGTSVL AQANQVPQNV LSLLRLEHHH HHH
Endotoxin	Less than 1 EU/g of rFlagellin, His as determined by LAL method.
Purity	≥ 95 % by SDS-PAGE analysis.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. - 12 months from date of receipt, -20 to -70 °C as supplied. - 1 month, 2 to 8 °C under sterile conditions after reconstitution. - 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at ≤ -20 °C. Further dilutions should be made in appropriate buffered solutions.