

mAb Standard, Mass Spec Grade

Cat. No. mAb Standard-57 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Description	Secukinumab monoclonal antibody is an interleukin-17A (IL-17A) inhibitor, which can be used for intact protein quality control (QC) and charge variants (CV). Also, it assists with the mAb subunit and N-glycan LC-MS analysis, purification, and system suitability verification of experimental methods.
Form	Mass spectrometry modified trypsin is supplied by lyophilized powder.
Molecular Mass	147940.0 Da
AASequence	Secukinumab Heavy Chain (CAS 875356-43-7): EVQLVESGGGLVQPGGSLRLSCA ASGFTFSNYWMNWVRQAPGKGLEWVAAINQDGSEKYYVGSVKGRFTISRDNKNS LYLQMNSLRVEDTAVYYCVRDYYDILTDYYIHYYFDLWGRGLTVTVSSASTKGPSV FPLAPSSKSTSGGTAALGCLVKDYFPEPVTVSWNSGALTSGVHTFPAVLQSSGLY SLSSVTVTPSSSLGTQTYICNVNHKPSNTKVDKRVKPKSCDKTHTCPPCPAPELLGG PSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPRE EQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYT LPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFL YSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK Secukinumab Light Chain (CAS 875356-44-8): EIVLTQSPGTLSPGERATLSCRASQSVSSSYLAWYQQ KPGQAPRLLIYGASSRATGIPDRFSGSGSGTDFTLTISRLEPEDFAVYYCQQYGSSP CTFGQGRLEIKRTVAAPSVFIFPPSDEQLKSGTASVVCLLNNFYPREAKVQWKVDN ALQSGNSQESVTEQDSKSTYLSSTLTLSKADYEKHKVYACEVTHQGLSSPVTKS FNRGEC

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Purity	99.5 % peak area measured by HPLC at 280 nm.
Storage	Store the lyophilized powder at -20 centigrade.
Shelf life	24 months at -20 centigrade.
Usage Notes	CV and Subunit by HPLC: The mAb antibody samples were digested with IdeS or IdeZ. Then, the consistency was analyzed by HPLC with SCX 2.1*50 mm column at 30 centigrade, followed by elution with mAb PH ion mobility buffer. N-Glycans by HPLC: 15 µg sample was reduced at 95 centigrade, denatured with TCEP, and digested with rPNGase F glycosidase. N-glycans were derivatized into glycosamines with N-glycan labeling reagents. Samples were analyzed by UPLC fluorescence detection using the HILIC column (FLR: Ex: 265 nm; Em: 425 nm).
Resuspension Buffer	Add a certain amount of pH Buffer to dissolve according to the experimental needs.