

293T whole cell lysate

Cat. No. 293T-2104H **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Species	Human
Tissue Type	Kidney
Preparation method	Culture cells in 6-well plate to 80% confluent, remove media and then rinse with cold 1x PBS. Scrape cells off the plate surface and pipette cells into microcentrifuge tube to spin down the cells. Remove PBS and then lysis and resuspend cells with 205 ul 1X sample buffer. Mix well and heat at 90°C for 10 minutes, and then cooled rapidly.
Size	100 uL
Recommended Usage	The lysate is ready to load on SDS-PAGE for Western blotting.
Storage Buffer	In ready-to-use 1X Sample Buffer (0.2 M Tris-HCl, 2% SDS, 10% glycerol, 0.2% 2-mercaptoethanol, 0.2% Bromophenol blue).
Applications	Western Blot;

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA