

Cibacron Blue Beads

Cat. No. Cibacron-25A **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Cibacron Blue Beads are primarily used for quick removal of albumins from protein-containing samples. The Cibacron Blue F3GA dye used in this product acts as a versatile binding agent binding to both charged and non-charged site in various molecules. Product was optimized for the removal of Bovine albumins while minimizing loss of other proteins present in sample.
Cat.no	Cibacron-25A
Matrix	Highly cross-linked spherical agarose, 4%
Functional group	Cibacron
Average particle size	45-165 µm
Dynamic binding capacity	0.8-1.1 mg/ml of gel
Recommended flow rate	1-3 ml/min/ cm ²

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA