

Heparin Beads

Cat. No. Heparin-27A **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Heparin Beads are made by coupling 20-35 KD size Heparin extracted from Porcine intestinal mucosa (Type I) to Sepharose beads. These beads are useful for isolating heparin binding protein such as Antithrombin III, lipoproteins, various enzymes and DNA binding proteins.
Cat.no	Heparin-27A
Matrix	Highly cross-linked spherical agarose, 4%
Functional group	Heparin
Average particle size	40-165 µm
Ligand	Heparin
Ligand concentration	0.7-1.2 mg/ml of 20-35 KD Heparin
Dynamic binding capacity	Thrombin binding capacity minimum 0.7 mg/ml of beads
Recommended flow rate	1-3 ml/min/ cm ²
pH working range	3-11

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA