

Native Bovine Collagen Type II

Cat. No. Collagen-56B **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Triplehelical domain of bovine nasal septum collagen type II. Gnd-HCl extraction, pepsin digestion, saltprecipitation, ion exchange chromatography, gelfiltration.
Species	Bovine
Source	Bovine Nasal Septum
Form	white, salt-free, strong hygroscopic. Lyophilized (0.5M HAC)
Quality Control Test	SDS and native polyacrylamide gel electrophoresis, collagenase digestion, Immunoblotting with anti-collagen, type I/II/III, rat (Rabbit).
Storage	Freezer
Reconstitution	Should be solubilized slowly in 10 mM acidic acid (0.2% w/v) and adjusted to neutral pH by adding 2/5 vol. 2.5x Carbonat buffer (0.25 M NaHCO ₃ , 1.25 M NaCl, 0.1% NaN ₃). Avoid freeze and thaw cycles. Fibril-forming molecules with the molecular formular: $\alpha 1(\text{II})_3$.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA