

DEAE polystyrene resin

Cat. No. diethylaminoethyl-006P **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Creative Biomart provides series of ion exchange chromatography media for high efficient separation and purification of a wide range of biomolecules such as proteins, oligonucleotides, and peptides. The ion exchange resins are composed of rigid, spherical, highly cross-linked polystyrene/divinylbenzene or polyacrylic beads.
Cat.no	diethylaminoethyl-006P
Matrix	cross-linked polystyrene
Functional group	diethylaminoethyl group: $-O-CH_2CH_2-N^+(C_2H_5)_2H$
Protein	diethylaminoethyl
Average particle size	35 μm
Dynamic binding capacity	20-60 mg/mL gel
Recommended flow rate	50-300 cm/h
Viscosity	The resin should be stored at 4-35 °C in 20 % ethanol
Chemical stability	Good chemical stability allows the use of the prepacked columns with most commonly used buffers and organic solvents, like for example 0.5 M NaOH, 50% ethanol, 20% isopropanol, 6 M GuHCl and 8 M urea.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Physical stability	Negligible volume variation due to changes in pH or ionic strength
pH working range	2-13
pH CIP range	1-14

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA