

Mini Filter Funnels for 10 mm and 19 mm flat widths, 5 pc

Cat. No. CBMart-24 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Mini Filter Funnels for 10 mm and 19 mm flat widths, 5 pc.

Description

Make loading and unloading dialysis tubing easy: Cut tubular membrane to length and clamp at one end, slip the open end onto the polypropylene funnel and secure with the silicone loop fastener. Open snap cap, fill with sample, and hang in container. Funnels are supported by a three-notched universal hanger that fits over the edge of common vessels.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA