

Inclusion Body Protein Isolation Maxi Kit

Cat. No. Inclusion Body-37A **Lot. No.** (See product label)

SPECIFICATION

Product Overview

The Inclusion Body Isolation Maxi Kit is a complete kit that is designed for the purification of recombinant proteins expressed in bacteria as inclusion bodies. Bacteria are widely used for the expression of different recombinant protein, however 70-80% of the proteins expressed in E. coli by recombinant techniques are contained in inclusion bodies (i.e., protein aggregates). This kit is stream-lined for rapid and efficient cell lysis, protein solubilization and on-column protein purification of up to 12 mg of proteins from 100 mL of bacterial culture.

Cat.no

Inclusion Body-37A

Description

The Inclusion Body Isolation Maxi Kit uses a proprietary Cell Lysis Reagent to selectively lyse the cells and release inclusion bodies in their solid form. Next, inclusion bodies are dissolved and their contents released using the IB Solubilization Reagent. Inclusion body proteins are then further purified by loading onto spin columns containing Norgen's proprietary resin (BIND). Non-specifically bound components in the solution can then be washed from the column, leaving the inclusion body protein bound to the resin (WASH). These specific proteins can then be eluted into the supplied Elution Buffer (ELUTE). Please see the procedure flowchart to the right.

Storage

The reagent should be stored at 4°C upon receipt of the kit. All other unopened solutions should be stored at room temperature. Once opened, the solutions should be stored at 4°C when not in use, except for the Basic and Acidic Binding Buffers, which should be stored at room temperature.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA