

Calcein AM Cell Viability Assay Kit (Fluorometric)

Cat. No. Kit-1114 Lot. No. (See product label)

SPECIFICATION

Product Overview

Description

Quantification of number of viable cells is an indispensable tool in Cell Biology research. Calcein AM Cell Viability Assay Kit is a fluorometric method for extremely sensitive quantification of viable cells that can detect as low as 50 viable cells in less than 30 min. Calcein AM is a non-fluorescent, hydrophobic compound that easily penetrates intact and; Live cells. Hydrolysis of Calcein AM by intracellular esterase produces a hydrophilic, strongly fluorescent compound that is retained in the cell cytoplasm and can be measured at Ex/Em = 485/530 nm. The measured fluorescence intensity is proportional to the number of viable cells. This assay kit provides an easy-to-use, non-radioactive, and high-throughput method for cell proliferation, cell viability, chemotaxis, cytotoxicity and apoptosis.

Applications

Measurement of cell proliferation in response to growth factors, cytokines, mitogens and nutrients

Analysis of cytotoxic and cytostatic compounds such as anticancer drugs, toxic agents and other pharmaceuticals

Assessment of physiological mediators and antibodies that affect cell growth.

Storage -20°C

Shipping Gel Pack

Size 1000 assays

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Kit Components	Calcein AM; Calcein Dilution Buffer
Detection method	Fluorescence (Ex/Em = 485/530 nm)
Features & Benefits	Convenient & non-radioactive; High-throughput; Highly-sensitive: detect as low as 50 viable cells in less than 30 min. Suitable for cell proliferation, cell viability, chemotaxis, cytotoxicity and apoptosis

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA