

Carbonic Anhydrase (CA) Activity Assay Kit (Colorimetric)

Cat. No. Kit-1116 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Description

Carbonic anhydrases (CA) (CA, 4.2.1.1) are zinc enzymes present in both prokaryotes and eukaryotes. They catalyze efficiently the reversible hydration of CO₂ to bicarbonate. Their important patho-physiological roles in respiration, pH and CO₂ homeostasis, secretion, gluconeogenesis, ureagenesis etc. make it an important drug target. Characteristic increase and decrease of CA Activity is observed in different physiological conditions and diseases such as anemia, Thalassemia, Hypothyroidism, Hyperthyroidism, in several cases of lung and liver diseases and Leukemia. CA Activity Assay Kit can be used to measure CA activity in biological samples like serum and hemolysates. The assay utilizes the esterase activity of an active CA on an ester substrate which releases a chromophore. The released product can be easily quantified using an absorbance microplate reader. In the presence of a CA specific inhibitor, the enzyme loses its activity which results in a decrease of absorbance. This assay kit is simple and can be used to measure CA activity in biological samples like hemolysate, serum and purified CA activity, CA activity in biological samples like hemolysate and serum in a high-throughput format. To measure the Specific CA activity we have used the specific inhibitor Acetazolamide, which is a potent inhibitor of CA I, CA II, CA IV, CA IX, CA XII etc. except CA III.

Applications

Detection of CA activity in hemolysates, serum;
Determination of enzymatic activity of purified CA

Storage

-20°C

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Shipping	Gel Pack
Size	100 assays
Kit Components	CA Assay Buffer; CA Dilution Probe; CA Positive Control; CA Substrate; CA Inhibitor (20 mM Acetazolamide); Nitrophenol Standard (2 mM)
Target Species	Mammalian
Detection method	Detection method - Colorimetric (OD 405 nm)
Features & Benefits	Simple, highly sensitive, high-throughput compatible; Easy protocol for the detection and quantification of active CA; Kit includes Acetazolamide, a CA inhibitor, and a stable, recombinant human CA

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA