

Luminescent ATP Detection Assay Kit

Cat. No. Kit-0116 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Luminescent ATP Detection Assay Kit irreversibly inactivates ATP degrading enzymes (ATPases) during the lysis step, ensuring that the luminescent signal obtained truly corresponds to the endogenous levels of ATP.
Size	300 tests
Description	Total levels of cellular ATP can be used to assess cell viability, cell proliferation and cytotoxicity of a wide range of compounds and biological response modifiers.
Applications	Functional Studies
Storage	Store kit at 4°C in the dark immediately upon receipt. Kit has a storage time of 1 year from receipt, providing components have not been reconstituted. Refer to list of materials supplied for storage conditions of individual components. Observe the storage conditions for individual prepared components in the Materials Supplied section. Aliquot components in working volumes before storing at the recommended temperature.
Kit Components	Detergent: 1 x 20ml; 4°C (Before Preparation); 4°C (After Preparation) Lyophilized ATP standard: 1 vial; 4°C (Before Preparation); 4°C (After Preparation) Lyophilized substrate: 3 vials; 4°C (Before Preparation); -20°C (After Preparation) Substrate Buffer: 1 x 20ml; 4°C (Before Preparation); -20°C (After Preparation)
Detection method	Luminescent

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Features & Benefits

Long luminescence signal: The half life of the luminescent signal is greater than 5 hours and therefore a special luminometer with injectors is not required. Fast: Results are obtained in less than 30 minutes. Simple: There are no separation steps and there are only 2 reagent additions. Homogenous: No cell harvesting or centrifugation required. Sensitive: The limit of detection is 5 cells in 100 μ . Wide linear dynamic range: Approximately from 1pM to 1M of ATP (this will vary depending on instrument sensitivity).

Assay time

0h 30m

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA