

Caspase-10 Fluorometric Assay Kit

Cat. No. Kit-0157 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Fluorometric assay kit for caspase-10 activity
Size	100 assay
Description	Activation of ICE-family proteases/caspases initiates apoptosis in mammalian cells. The Caspase-10 Fluorometric Assay Kit provides a simple and convenient means for assaying the activity of caspases that recognize the sequence AEVD. The assay is based on detection of cleavage of substrate AEVD-AFC (AFC: 7-amino-4-trifluoromethyl coumarin). AEVD-AFC emits blue light ($\lambda_{\text{max}} = 400 \text{ nm}$); upon cleavage of the substrate by caspase-10 or related caspases, free AFC emits a yellow-green fluorescence ($\lambda_{\text{max}} = 505 \text{ nm}$), which can be quantified using a fluorometer or a fluorescence microtiter plate reader.
Applications	Detect early/middle stages of apoptosis; differentiate apoptosis from necrosis.
Target Species	Mammal
Usage	For research use only. Not intended for any diagnostic or therapeutic purpose.
Storage	Store at -20°C . Product has an expected shelf life of 6-12 months.
Kit Components	Cell Lysis Buffer; 2X Reaction Buffer; AEVD-AFC (1 mM); DTT (1 M)
Detection method	Fluorescence microtiter plate reader (400-nm excitation filter and 505-nm emission filter) and fluorescence microscopy

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

**Compatible Sample
Types**

Cell and tissue lysates

Features & Benefits

1. Simple one-step procedure; takes only 1-2 hours. 2. Fast and convenient. 3. Comparison of the fluorescence of AFC from an apoptotic sample with an uninduced control allows determination of the fold increase in AEVD-dependent caspase activity.

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA